


LEADING EDGE OF DR SOLUTIONS

myvet DR Precision



ENHANCE YOUR IMAGING EXPERIENCE WITH LEADING EDGE INTEGRATION

- Cutting Edge of Imaging Software Technology
- DR Integration with Automatic X-Ray Exposure Detection
- Workstation with Fully Integration to The Existing System
- Auto Triggering Technology



Designed to offer new levels of digital radiographic image quality without the need to modify existing x-ray equipment

◀ Digital Radiography Flat Panel Detector

MyVet DR 1717 features auto-trigger signal sensing technology. Within seconds of image preview time, patient workflow and overall productivity are improved and patient latency is reduced, providing an economical and ideal solution for x-ray diagnostics.

1717 Tethered FPD

- Fixed cassette sized: No modifications required
- Image quality: 127 μm , Cesium-Iodide
- 1GB connection speed
- Auto triggering technology



◀ Acquisition Workstation and software

VetView Pro is a user-friendly veterinary software that streamlines image acquisition, processing, and management in real time

- DICOM compliant System
- Premier Image Processing with HDR
- Records Dose History
- Measure the Length and Angles of the Image
- Add annotation text, graphics and markers to an image



◀ Smart Viewer with Diagnostic Tool

VetView Pro features an intuitive interface with customizable tools, enabling veterinarians to easily adjust images and tailor the software to their specific needs.

